# HITTEPIT<sup>®</sup>

Treets

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# The Natural Cherry Pit Heating Pad

Hittepit is a unique and natural product that can be used both hot and cold.

Hittepit can be heated up in the microwave. It retains the warmth for a long time and releases it very slowly to the surrounding area. In contrast to a hot water bottle, Hittepit does not draw heat from the body after it has cooled down. Hittepit conforms to the shape of the body and is a 100% natural product.

Use Hittepit where heat brings relief, for muscle or stomach ache, back complaints or menstruation pains.

Hittepit can also be used to help relax a stiff neck for example, and for a hand or foot massage and where muscle movement requires activating.

**Recommendation:** you can also use your cherry pit cushion as a cold compress. Put a plastic bag around the cushion before putting it in the freezer compartment or the freezer. The cushion will provide ideal cooling after a few hours.

## Instructions for use

The cherry pit cushion should be heated up in the microwave for a maximum of 2.5 minutes at 700 W. Adjust the heating time if your microwave uses a different power level.

## Power in Watts

700 W

800 W

900 W

1,000 W

## Heating time

Max. 2.5 min.

Max. 2 min.

Max. 1.5 min.

Max. 1 min.

## 12 steps for using your cherry pit cushion in the right way

- 1) Make sure that the microwave is clean and fat free before heating up the cherry pit cushion.
- 2) Check the power level of the microwave.
- 3) Set the heating time according to the power level of the microwave as shown in the table above.
- 4) Remove all covers or towels from around the cushion.
- 5) Moisten the cherry pit cushion with water. It is better for the cushion to be too wet than too dry!

- 6) Put the cushion on a plate in the middle of the microwave. Make sure that the cushion is not touching the sides of the microwave and that the turntable is not blocked.
- 7) Make sure that no other functions, such as the grill or hot air, are switched on at the same time. Always remove any grating from the microwave before heating up the Hittepit
- 8) Stay near to the microwave when you are heating the cushion and keep an eye on the heating time.
- 9) Turn the cushion over halfway through the heating time.
- 10) Check the temperature of the cushion before using it, both when used for heating and cooling. If the cushion does not feel hot when first touched shake it for a little while.
- 11) Never reheat the cushion immediately after heating it.
- 12) Wait at least two hours before heating the cushion again and check that it has cooled off properly before doing this.

## Important points

Read the following instructions carefully before using the cherry pit cushion and make sure that you adhere to these to guarantee safety and a great deal of pleasure in using the Hittepit.

Following the instructions will prevent fire and other damage due to overheating and drying out.

- This product is only suitable for the microwave
- Do not allow children or persons without experience to heat up the cushion.
- Select the right power and warming up time in the microwave carefully
- Keep to the instructions and the 12 steps for use
- Use the cherry pit cushion to warm up a bed or cot but do not leave the cushion in a bed or cot with a baby or child in it
- Prevent the cherry pit cushion from drying out by moistening it with water each time before use. It is better for the cushion to be too wet than too dry!

- The working of the microwave spreads the heat evenly from the inside to the outside. Because of this the cherry pit cushion does not feel warm to the touch when you first take it out of the microwave. Shake the cushion and you will notice that it really is warm.
- Every microwave can have hotspots. Turn the cushion around halfway through the heating time to ensure the best possible distribution of energy and even heating
- Allow the cushion to cool down for two hours before it can be heated up again.
- Do not use the cherry pit cushion any more if it has been damaged by overheating or other misuse
- Exceeding the warming up time or using a power level that is too high (Watts) will not result in heat being given out for longer but can lead to damage by scorching or burning
- This cushion must never be used unsupervised by and for babies and young children or persons with limited reactions
- Do not use if there are symptoms of inflammation. Persons with diabetes, blood circulation, sensorial or neurological problems or people who are sensitive to temperatures should not use this product
- Do not use Hittepit as a substitute for medical treatment
- Always keep these instructions for use together with your cherry pit cushion
- Do not wash the cushion in the washing machine  
If the Hittepit has become dirty, for example by fat, then it is damaged and should no longer be used. Throw the cushion away immediately!

Make sure you consistently follow the times and Wattages listed above and always check the Hittepit's temperature before using it. Nature's Choice B.V. is not responsible for any damage resulting from improper use or failure to follow the Hittepit user instructions. Keep these instructions near your cushion!

**NB: IF THE INSTRUCTIONS FOR USE ARE NOT FOLLOWED THE PRODUCT CAN CATCH FIRE**

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Copyright © Nature's choice B.V. 2009. Alle rechten voorbehouden.